

DermaLume 2X[®]

Spot and Scalp Treatment

OVERVIEW

The DermaLume 2X[™] by National Biological is a handheld phototherapy light wand. It delivers localized treatment to small areas of psoriasis and vitiligo. For patients with scalp psoriasis, it comes with a distance comb. This comb has two important functions: 1) It enables the patient to separate hairs, allowing light to penetrate directly to the scalp and 2) It helps the patient maintain a safe, consistent distance between the device and his skin.

KEY FEATURES

- Ideal for scalp Psoriasis treatment or small patches of Psoriasis or Vitiligo
- Multiple wand positions for comfortable hold
- Safety features include automatic lamp shutoff at end of treatment
- Ships fully assembled
- Comes with carrying case for travel or storage and goggles
- Covered by many insurance plans

ELECTRICAL RATING

- 110-120 VAC, 60 Hz
- 0.4 A with IEC grounded plug

DIMENSIONS

- Output window: 3.3" x 5.4" / 17.5 sq in (8 cm x 14 cm / 113 cm²)
- Wand (folded): 4.1" x 7.5" x 2.8" (11 cm x 19 cm x 7 cm)
 - Wand (open): 4.1" w x 15" x 2.8" (11 cm x 38 cm x 7 cm)



LAMPS

- Two 311nm Narrowband UVB lamps, 18 watts UV energy

INTERNATIONAL PRODUCT SPECS

- For countries with electrical current other than U.S. standard, DermaLume 2X[™] device will require a step-down transformer.

WARRANTY

- All parts (excluding lamps) and labor carry a 1-year full replacement warranty
- Lamps carry a 90-day warranty



Authorized Distributor

Dermitech Phototherapy
(214) 377-8144
sales@dermitech.com