

Panosol 3D® Full Body Treatment



OVERVIEW

The Panosol 3D® by National Biological provides full-body treatment in the privacy of your home. A single panel with two adjustable wings, the Panosol 3D enables the patient to treat the entire body in two sessions instead of four by allowing exposure of the sides of the body to UV light at the same time as the front or back. Wings close for easy storage.

KEY FEATURES

- Multi-panel design reduces treatment time by treating three sides of the body at once
- High-efficiency reflector design maximizes lamp output
- Wings close for easy storage
- Operates on standard household current
- Covered by insurance plans

ELECTRICAL RATING

- 110-120 VAC, 50 /60 Hz, grounded plug
- NB UVB device-8 A

DIMENSIONS (INCLUDING FEET)

- Closed: 29.5" W x 81" H x 15" D – 73 cm x 206 cm x 38 cm
- Open as used: 48.5" W x 81" H x 21" D – 123 cm x 206 cm x 53 cm

LAMP OPTIONS

- 10 Narrowband UVB 311nm lamps
- 10 UVA lamps
- 10 UVA1 lamps

INTERNATIONAL PRODUCT SPECS

- For countries with electrical current other than U.S. standard, the Panosol 3D will require a step-down transformer.

WARRANTY

- All parts (excluding lamps) and labor include a 3-year full replacement warranty
- Lamps carry a 90-day warranty



Authorized Distributor

Dermitech Phototherapy
(214) 377-8144
sales@dermitech.com