

# Panosol II® Full Body Panel



## OVERVIEW

The Panosol II® by National Biological is a series of light panels that come in two heights: 2 feet and 6 feet.

## KEY FEATURES

- High-efficiency reflectors maximize lamp output
- Safety features include automatic lamp shutoff at end of treatment
- 6' model comes with casters
- Operates on standard house current
- Comes with goggles
- Covered by many insurance plans

## ELECTRICAL RATING

- 110-120 VAC, 50 /60 Hz, grounded plug
- 5A maximum current

## DIMENSIONS (INCLUDING FEET)

- 6-foot: 29.5" W x 15" D x 78" H— 75 cm x 38 cm x 198 cm
- 2-foot: 29.5" W x 11" D x 29.5" H— 75 cm x 28 cm x 75 cm

## LAMP OPTIONS

- 6-foot:
  - 6 Narrowband UVB 311nm lamps
  - 6 UVA lamps
- 2-foot
  - 8 Narrowband UVB 311nm lamps
  - 8 UVA lamps

## INTERNATIONAL PRODUCT SPECS

- For countries with electrical current other than U.S. standard, the Panosol II will require a step-down transformer.

## WARRANTY

- All parts (excluding lamps) and labor include a 3-year full replacement warranty
- Lamps carry a 90-day warranty



**Authorized Distributor**

**Dermitech Phototherapy**  
**(214) 377-8144**  
**sales@dermitech.com**